# CALBRATORE PER CHIMICA CLINICA / CLINICAL CHEMISTRY CALIBRATOR

Calibratore a base umana / Human based calibrator.

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### BIBLIOGRAFIA

- Department of Labor, Occupational Safety and Health Administra 29 CFR Part 1910.1030 Fed. Register Juli 1, 1998; 6:267-280
- Council Directive 90/679/EEC. Official journal of the European directives N°.L374 from Dec.31,1990:1-12

## PRODUCT CHARACTERISTICS

Is a human multiparametric calibration serum. The concentrations and the activities have been selected to ensure optimum calibration of automatic analyzers. Particularly, these values were determined in our laboratory on COBAS MIRA S

### REAGENTS

Human serum with chemical additives and tissue extract of human and animal origin. The origin of biological additives is as follow:

AST/GOT, ALT/GPT, LDH:	porcine heart
ALBUMIN, CHOLESTEROL:	bovine plasma
ALDOLASE, CREATINE KINASE:	rabbit muscle
AMYLASE, TOTAL, PANCREATIC; LIPASE:	porcine pancreas
CHOLINESTERASE:	human serum
ACID PHOSPHATASE:	human prostate / potato
ALKALINE PHOSPHATASE:	placenta (human, recombinant)
GAMMA-GT:	porcine kidney
GLDH:	bovine liver
TRANSFERRIN	human Cohn IV fraction
TRIGLYCERIDES	chicken egg yolk

The concentrations and the activities are lot specific. The exact values are in the enclosed values sheet.

### PRECAUTIONS AND WARNING

For in vitro diagnostic use only.

*The product is obtained using only blood of donors tested by an FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However as no test method can rule out the potential risk of infection with absolute certainty, the material should be handled as potentially infectious, just as carefully as a patient sample.*

### WASTE DISPOSAL

Please consult local regulations for a correct waste disposal.

S56: Dispose of this material and its container at hazardous or special waste collection point.

S57: Use appropriate container to avoid environmental contamination.

S61: Avoid release to the environment. Refer to special instructions / Safety data sheets.

### PREPARATION

Carefully open one bottle of calibrator, avoiding the loss of lyophilized, and pipette exactly 3.0 ml of distilled water. Carefully close the bottle and dissolve the contents completely by occasional gentle swirling within 15 minutes.

Avoid the formation of foam.

### STORAGE AND STABILITY

Store at +2-8°C. Lyophilized calibrator serum is stable at +2-8°C until the stated expiration date.

Stability of the components in the reconstituted calibrator:

at 15-25°C	8 hours
a 2-8°C	2 days
a (-15) - (-25)°C	2 weeks (when frozen once)

Stability of Bilirubin in the reconstituted calibrator (when stored protected from light):

at 15-25°C	4 hours
a 2-8°C	1 day
a (-15) - (-25)°C	2 weeks (when frozen once)

Stability of Acid Phosphatase in the reconstituted calibrator (protected from light):

at 15-25°C	4 hours
a 2-8°C	1 day
a (-15) - (-25)°C	2 weeks (when frozen once)

Store reconstituted calibrator tightly capped when not in use.

### REQUIRED MATERIALS NOT PROVIDED

Automatic micropipette, distilled water, general laboratory equipment.

### ASSAY

Dispense the required volume into a sample cup and use as advice.

### RESULTS

The calibration values were determined using the method mentioned in the enclosed value sheet. Determinations were performed under strictly standardized conditions on Roche / Hitachi analyzers using Roche system reagents and Roche master calibrators. The calibration values were obtained via single assays performed in different laboratories in several independent series. The calibration value specified is the median of the values obtained. Giesse Diagnostics calibration takes place in own laboratory on COBAS MIRA S analyzer and with own reagents.

### REFERENCES

- Department of Labor, Occupational Safety and Health Administra 29 CFR Part 1910.1030 Fed. Register Juli 1, 1998; 6:267-280
- Council Directive 90/679/EEC. Official journal of the European directives N°.L374 from Dec.31,1990:1-12